IVD

REF Q-1 Pro

qLabs[®] ElectroMeter

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Tel.: +86 755 86296766

**Micropoint Biotechnologies, Inc.**

6F, No.3, Industry 5 Road, Shekou,
Shenzhen, Guangdong 518067, China

Tel.: +86 755 86296766

Fax: +86 755 86673903

www.micropointbio.com

**OBELIS S.A**

Bd. General Wahis 53,
1030 Brüssel, Belgien

Tel.: +32 2 732 59 54

Fax: +32 2 732 60 03

Table of Contents

1. Introduction.....	3
1.1 Before You Start	
1.2 Intended Use	
1.3 Test Principle	
1.4 Packaging	
2. Accessory List	4
3. Product Overview	4
4. Operation Summary	6
5. Before Testing	6
6. Precautions, Limitations and Warnings	7
6.1 Care of Your qLabs® ElectroMeter	
6.2 Patient Health Status	
6.3 Performing a Test	
6.4 Collecting a Fingerstick Blood Sample	
6.5 Electromagnetic Compatibility	
7. Powering the qLabs® ElectroMeter ON and OFF	9
7.1 Power Connection	
7.2 Powering the qLabs® ElectroMeter On	
7.3 Powering the qLabs® ElectroMeter Off	
8. qLabs® ElectroMeter Settings	11
8.1 Enter Setup Mode	
8.2 Setting the date	
8.3 Setting the Time	
8.4 Entering Patient ID	
8.5 Bluetooth Settings	
8.6 Lot Number Setup	
8.7 Setting a Target Range For Your INR Results	
8.8 Selecting Displayed Units	
8.9 Selecting Language	
9. PT Test	16
9.1 Entering Test Mode	
9.2 Insert a Test Strip	

9.3 Entering the Strip Code and Lot Number (if enable)	
9.4 After confirming the Code	
9.5 Heating Up	
9.6 Getting a Fingerstick Sample	
9.7 Performing the Test	
9.8 Results and Explanation	
9.9 Finishing the Test	
10. Reviewing Previous Results	22
10.1 Reviewing the Previous Results	
10.2 Fast Deleting the Previous Results	
11. Quality Control	26
11.1 Onboard QC Test	
11.2 Liquid QC Test	
12. Maintenance	27
12.1 Care and Cleaning of Your Screen	
12.2 Installing/Replacing Batteries	
12.3 Servicing	
13. Troubleshooting	28
14. Symbols	30
15. Performance Characteristics & Product Specifications	31
15.1 Operating Condition	
15.2 Product Specifications	
16. Warranty	32
17. Contact Information	34

1. Introduction

1.1 Before You Start

Before using the qLabs® ElectroMeter to test Prothrombin Time (PT) and International Normalized Ratio (INR), read the entire User's Manual carefully.

1.2 Intended Use

The qLabs® PT-INR monitoring system is used for quantitative measurement of Prothrombin Time in fresh, capillary whole blood samples. The qLabs® PT-INR monitoring system is intended for in vitro diagnostics and is not intended to be used for screening purposes.

1.3 Test Principle

The qLabs® PT-INR monitoring system consists of the qLabs® ElectroMeter and Prothrombin Time (PT)-International Normalized Ratio (INR) test strips. The qLabs® ElectroMeter automatically detects the insertion of a qLabs® PT-INR test strip and heats the strip to a preset operating temperature. When a drop of blood is added to the sample well of the strip, it flows through channels to two reaction zones: a Test Zone and a Control Zone. In these zones the blood mixes with pre-printed reagents and begins to coagulate. Each reaction zone contains a pair of electrodes to which a constant voltage is applied by the qLabs® ElectroMeter. As the coagulation of the blood proceeds, the current monitored across the two electrodes changes. The qLabs® ElectroMeter detects the change of the current and determines the PT and INR results.

1.4 Packaging

The qLabs® ElectroMeter is packed individually.

Upon receiving the package, please open and remove the packaging materials. Place the qLabs® ElectroMeter on a flat surface and connect the power adapter.

Do not keep the qLabs® ElectroMeter in direct sunlight, near a high heat source, or near an area with a strong magnetic field.

2. Accessory List

The qLabs® ElectroMeter includes the following accessories:


- qLabs® ElectroMeter
- Power Adapter
- User's Manual
- Quick Start Guide
- Lancet Device
- Lancet Needles
- Carrying Case
- Instructional DVD
- qLabs® Data Manager Installation DVD
- qLabs® Data Manager Installation Guide
- qLabs® eCable

NOTES:

- To avoid the possibility of damage, use only the power adapter supplied with the qLabs® ElectroMeter .
- Standard AA batteries can also be used for the qLabs® ElectroMeter.
- The qLabs® eStation and the qLabs® Data Manager can be purchased separately from Micropoint Biotechnologies, Inc. Please contact your local distributor for ordering information.

3. Product Overview

Identification of Components



Front View

1. LCD Screen
2. Return/Power Button
3. Right Button
4. Left Button
5. OK Button
6. Test Strip Guide

Back View

1. Battery Cover
2. Battery Cover Screw
3. Battery Cover Release


Top View


1. Power Supply Jack
2. Data Port

qLabs® PT-INR Test Strip

1. Sample Well
2. Test Zone
3. Electrode Leads

4. Operation Summary

5. Before Testing

Before using the qLabs® monitoring system to do a test, you will need:

- qLabs® PT-INR test strip
- Power supply or batteries
- Alcohol wipes
- Gauze and Band-Aid
- Lancet device
- Puncture proof container

6. Precautions, Limitations and Warnings

The qLabs® PT-INR monitoring system is intended for in vitro diagnostics use only. Before using this system to test PT and INR, take special note of CAUTIONS throughout this User's Manual.

6.1 Care of Your qLabs® ElectroMeter

- The qLabs® ElectroMeter is a delicate instrument and should be handled with care. Dropping or other mishandling may cause malfunction of the qLabs® ElectroMeter.
- The qLabs® ElectroMeter should be transported in a carrying case or a secure container.
- DO NOT spill any liquid on the qLabs® ElectroMeter. If this should occur, immediately contact your local distributor from Micropoint Biotechnologies, Inc.
- DO NOT store the qLabs® ElectroMeter below -10 °C or above 40 °C.
- DO NOT use the qLabs® ElectroMeter for any other types of test strips not provided by Micropoint Biotechnologies, Inc.

6.2 Patient Health Status

Current patient health status may cause inaccurate or unexpected test results. It is important to take certain health factors into consideration when interpreting the test results and deciding on a course of action for your patients. Failure to do so may lead to an incorrect interpretation of the PT-INR result.

6.3 Performing a Test

- The qLabs® ElectroMeter should be operated on a level surface that is free of vibration. Testing on an uneven or unstable surface may cause inaccurate results. DO NOT hold the qLabs® ElectroMeter in your hands during the testing.
- The blood sample must be applied to the test strip immediately after collection. Otherwise, the blood sample may begin clotting, causing inaccurate results.

- Plasma samples or whole blood samples with anticoagulation reagents are not suitable for the qLabs® PT-INR test.
- DO NOT move or touch the qLabs® ElectroMeter during testing.

During Testing DO NOT switch the power supply from the adapter to the battery or the vice versa!

6.4 Collecting a Fingertick Blood Sample

- Sanitize the site of the fingertick with alcohol wipe before testing.
- Use a fresh capillary fingertick blood sample for testing.
- Squeezing the fingertick site excessively (milking) may release interstitial fluid into the blood sample, which may lead to inaccurate results.
- The site of the fingertick must be completely dry. If any alcohol remains on the finger, it may cause hemolysis, which may lead to inaccurate results.

6.5 Electromagnetic Compatibility

- This equipment has been tested and found to comply with applicable EMC emission requirements as specified in EN 61326-1: 2006 and EN 61326-2-6: 2006, immunity to electrostatic discharge as specified in IEC 61000-4-2, and immunity to radio-frequency interference at the frequency range and test levels specified in IEC 61000-4-3. The emissions of the energy used are low and not likely to cause interference in nearby electronic equipment.
- Do not use this instrument in a dry environment, especially if synthetic materials are present. Synthetic clothes, carpets, etc., may cause damaging static discharges in a dry environment. This may cause incorrect operation or damage to the device. The recommended humidity operating range for qLabs® ElectroMeter is 10% to 90%.
- Do not use this instrument near cellular or cordless telephones, walkie talkies, garage door openers, radio transmitters, or other electronic equipment that are sources of electromagnetic radiation, as these may interfere with the proper operation of the instrument.

7. Powering the qLabs® ElectroMeter ON and OFF


7.1 Power Connection

The qLabs® ElectroMeter comes with a power supply which must be plugged into the qLabs® ElectroMeter power supply jack. Alternatively, the qLabs® ElectroMeter can use 4 standard AA batteries (See Section 12.2).

The messages on the four-cell grid displayed at the bottom of the screen provide hints for the functions of the buttons.

7.2 Powering the qLabs® ElectroMeter On


When the qLabs® ElectroMeter is initially powered on, it will beep and prompt to insert a test strip automatically.

When the qLabs® ElectroMeter is in the standby mode, you have two choices to power on the monitor:

1. Press and hold the “Power” button for 2 seconds to turn on the qLabs® ElectroMeter, or
2. Insert a fresh test strip if you are ready to perform a test. The qLabs® ElectroMeter will then prompt you to confirm the strip code.

Turning the qLabs® ElectroMeter ON in CHOICE 1, allows the user to:

- Perform a test, or
- check the version of the software and the manufacturer information by pressing the “VER” button, or

- enter the Setup Mode by pressing the “SET” button to change the ElectroMeter setting, or
- enter the Memory mode by pressing the “MEM” button to review the previous results.


Turning the qLabs® ElectroMeter ON in CHOICE 2, allows the user to:

- Perform a test.


See the appropriate sections of this User’s Manual for complete information for all of these options.

7.3 Powering the qLabs® ElectroMeter Off

In the interface where the monitor prompts you to insert a test strip, press the functional “OFF” button to enter the shut down mode.

In shut down mode, the user can:

- Press the functional “OK” button to turn off and enter standby mode, or
- press functional “BACK” button to return to the previous interface.


NOTES:

- Manually turn off the qLabs® ElectroMeter if not in use to conserve the batteries.
- If the qLabs® ElectroMeter is left idle for about 5 minutes, it will turn off automatically and enter standby mode.


8. qLabs® ElectroMeter Settings

Follow the steps below to set up the qLabs® ElectroMeter.

8.1 Enter Setup Mode


When the qLabs® ElectroMeter is turned on for the first time, the display screen will prompt the user to insert a test strip automatically.


When the screen prompts the user to insert a test strip, press the functional “SET” button to enter the Setup mode.

8.2 Setting the date

The date format is DD/MM/YY. The first set of digit (YY) will be highlighted. Use the functional “+” and “-” buttons to change the date.

Press the functional “NEXT” button to go to the second digit (MM) after a correct date is selected. Change the month and the day using the same method.

8.3 Setting the Time

The time is in 24-hour format: "HH:MM" (Hour: Minute).

01/02/12		
12:00		
DD/MM/YY HH:MM		
-	NEXT	+
	OK	

After setting the day, the setting option will prompt you to the hour setting mode. The hour will be highlighted on the display. Use "+" and "-" functional buttons to adjust the hour. Press "NEXT" functional button when you have the correct hour on the display.

01/02/12		
12:00		
DD/MM/YY HH:MM		
-	NEXT	+
	OK	

The minute will be highlighted on the display. Use "+" and "-" buttons to change the minute. Press "NEXT" button when you have the correct minute on the display.

8.4 Entering Patient ID


PATIENT ID		
	NEXT	
	OK	

You can press "" to enter a new patient ID or change the saved patient ID. Press "NEXT" button to enter the next test, press "OK" to close the setting and return to the main menu. To distinguish the sources of data that uploaded into the centre database, the patient ID you input will be saved into every test record automatically.

PATIENT ID		
A000  _____		
-	NEXT	+
	OK	

The Patient ID can be composed of numbers or letters, with a length of no more than 16 digits. Use "+" and "-" buttons to change the highlighted digit. Press "NEXT" button to accept and advance to next digit. If the input of the Patient ID is completed, press "OK" button to accept.


8.5 Bluetooth Settings

Press “+” or “-” buttons to turn on/off the bluetooth. When the Bluetooth function is on, the user can use optional Bluetooth dongle to connect the meter to the user management software.

Please refer to the Data Manager manual for Bluetooth driver installation and usage.

8.6 Lot Number Setup

You can press the functional “+” and “-” buttons to enable or disable the LOT NUMBER input function. If the LOT NUMBER input function is enabled, you will be asked to enter the Lot Number before performing each test. The Lot Number contains the expiration date of the strip or liquid control. The meter will check whether the current date exceeds the expiration date of the strip or liquid control. If the current date exceeds the expiration date, the meter will report an error.

NOTE: Please make sure the date setting is correct if you want to use the input function of LOT NUMBER.

8.7 Setting a Target Range For Your INR Results

An INR target range may be set to alert the user when the INR results are out of range.

If results fall outside the target range, or if results fall within the target range but are not consistent with the patient's current state of health (e.g. they have symptoms such as bleeding or bruising) follow your usual procedures for taking corrective action.

Each physician should establish expected values for his or her patient population or individual patients. Differences in reagents, instruments and pre-analytical variables can affect prothrombin time results. These factors should be considered when comparing different prothrombin time test methods.

TARGET: ON 0.7 1.4		
-	NEXT OK	+

Use “+” and “-” buttons to switch between the Target ON and OFF modes. If “OFF” is selected and the “NEXT” button is pressed, the screen will prompt you to the date setting mode. If “ON” is selected and “NEXT” is pressed, refer to the following sections.

If you switch the Target ON,

TARGET: ON LO 0.7 1.4		
-	NEXT OK	+

press the “NEXT” button and the “LO” message will appear on the display, indicating the user to set the lower limit of the target range. Use the “+” and “-” buttons to change the lower limit. The “LO” message will also be flashing on your test result display whenever your result is below the lower limit.

Press the “NEXT” button when you have the correct lower INR limit and are ready to set the higher limit of the target range.

TARGET: ON 0.7 HI 1.4		
-	NEXT OK	+

The “HI” message will appear on the display. Use the “+” and “-” buttons to change the higher limit. The “HI” message will also be flashing on your test result display whenever your result is above the higher limit.

Press the “OK” button to accept any setting. The system will save your setting and will return to the display that prompts the user to insert a test strip.

If you switch the Target OFF,

TARGET: OFF 0.7 1.4		
-	NEXT OK	+

press the “OK” button. The system will save your settings and will return to the display that prompts the user to insert a test strip. The “LO” or “HI” message will not be flashing on your test result.

If you would like to reset the settings, press the “NEXT” button to go back to the date setting. Repeat section 8.2 until the desired setting is selected.

NOTES:

- In the Setup Mode, pressing the “OK” button will complete all setup steps and save your settings, then return to the display that prompts the user to insert a test strip.
- If the power is cut off, current settings will be lost and the qLabs® ElectroMeter will reset to the factory default settings. To save settings, install / replace the batteries while keeping the power adapter connected to the power supply.

8.8 Selecting Displayed Units

UNITS
 ✓ PT/INR/QC
 PT/INR
 INR

You can use the right button or the left button to select one of the three display modes, which are INR, PT/INR and PT/INR/QC.

8.9 Selecting Language

LANGUAGE
 ✓ ENGLISH
 FRANÇAIS
 DEUTSCH
 ITALIANO


You can use the right button “▶” or the left button “◀” to select the desired language and then press the “OK” button to return to the main menu.

NOTE:

- The language option may be different, depending on the software configuration.

9. PT Test

9.1 Entering Test Mode

When the qLabs® ElectroMeter is turned on by pressing the “OK” button, the display screen will prompt the user to insert a test strip. Go to Section 9.2 for test strip insertion.

When the qLabs® ElectroMeter is turned on by inserting a fresh test strip, the display screen will prompt the user to confirm the strip’s code. Go to Section 9.3 for entering the strip code.

9.2 Insert a Test Strip

Remove a fresh test strip from its foil pouch. Insert the strip into the test strip guide.

The following two observations indicate that you inserted the strip correctly:

1. The electrode leads go in first.
2. You can read “PT-INR” clearly on the blue area of the test strip from left to right.

9.3 Entering the Strip Code and Lot Number (if enable)

When the test strip is inserted into the qLabs® ElectroMeter, the display will prompt the user to enter a 7-digit Strip Code.

If the displayed code and lot number (if enable) match the code on the pouch, press the “OK” button to confirm. Otherwise, press “✎” button to edit the strip code and lot number (if enable).

STRIP CODE		
PT/INR		
1505320		
LOT: 3052E0503		
-	NEXT	+
	OK	

Use the “+” and “-” buttons to change the first digit that is highlighted.

STRIP CODE		
PT/INR		
1505320		
LOT: 3052E0503		
-	NEXT	+
	OK	

Press the “NEXT” button to accept and advance to the next digit. To change numbers more quickly, hold down the functional “+” and “-” buttons. When the 7-digit Strip Code is entered correctly, press the “OK” button to finish the editing and return to the code confirming interface.

NOTES:

- If the Lot Number input function is enabled, you will be asked to input the strip code and the lot number; otherwise, you just need to input the strip code.
- Always match the code on the display with the strip code on the strip pouch. Failure to do so may yield inaccurate results.

9.4 After confirming the Code


EXPIRED STRIP!		
	OK	

If you input the lot number, the meter will automatically check whether the current date exceeds the expiration date of the strip or liquid control after confirming the code is correctly entered. If exceeds, an error will be reported. As a result, you need to use strips within the expiration date to perform the test. If you didn't input the lot number, the meter will skip the expiration date inspection and start the test directly.

9.5 Heating Up

	HEATING
	

After confirming the code, the qLabs® ElectroMeter will heat up automatically for the test.


When qLabs® ElectroMeter is ready for sample testing, it will beep and prompt the user to apply a blood sample.

9.6 Getting a Fingerstick Sample

Good fingerstick technique is beneficial to both the professional user and the patient. Getting a good fingerstick will ensure minimal discomfort for the patient and accurate test results for the user. Please sanitize the site of the fingerstick with alcohol wipe before testing.


9.6.1 Maximize Good Blood Circulation

- 1) If needed, briefly warm the hand in warm water or with a heating pad.
- 2) Massage the finger with a downward motion several times before performing the fingerstick.
- 3) Lower the hand to below the heart level when collecting the blood drop.


9.6.2 Identify a site on the finger to puncture

- 1) On one of the middle fingers of either hand.
- 2) Near the top of the finger on either side.
- 3) Away from any calluses or scars.


9.6.3 Clean the selected area with 70% isopropyl alcohol, or an alcohol pad. Dry thoroughly with cotton or gauze.


9.6.4 Puncture the finger following the instructions for the lancet that you are using.

9.6.5 Apply gentle, continuous pressure until a large, hanging drop of blood forms.


9.7 Performing the Test


Apply the blood directly on the sample well of the strip. The minimum sample volume is 10 μ L.

After the blood sample is added, the qLabs® ElectroMeter will perform the test automatically.

Please add sufficient sample blood into the strip sample well, or the qLabs® ElectroMeter will not start the test and will remain in the "ADD BLOOD" screen. REMOVE the test strip if this screen stays on for more than 1 min. Go back to Section 9.2 and use a fresh test strip.

When testing is complete, the qLabs® ElectroMeter will beep and the test results (PT, QC and INR) will appear on the screen along with the date and time.

CODE: 1505320		
LOT: 3052E0503		
01/05/2014 14:30		
		
		

If the test results need to be reported instantly, press “→” button to print or upload the results.

		
		
SET	OFF VER	MEM

When the “↶” button is pressed, the qLabs® ElectroMeter will prompt the user to remove the test strip.

NOTES:

- Test results are automatically stored in the qLabs® ElectroMeter memory including all the information of the test (See Section 10 for reviewing past results).
- The qLabs® ElectroMeter memory can store up to 200 test results. If the memory is full, the new test result will overwrite the oldest test result.

9.8 Results and Explanation

• Normal Range:

Results for normal blood were determined by testing 120 subjects who were not taking anticoagulant medication. The ranges found were INR: 0.7-1.4.

• Therapeutic Range:

Therapeutic ranges are determined for each patient individually by their health care professional. While most recommendations are within an INR range of 2 to 4.5, values outside the indicated range may be encountered.

• Re-establish Therapeutic Target Range:

There could be a bias from one test system to another. You must re-establish the INR therapeutic target range at the initial use of qLabs® system. This range may be different from your previous therapeutic target range.

- **Unexpected Results:**

When the qLabs® ElectroMeter displays a PT-INR result outside of the expected therapeutic range, it may or may not be due to an unusual clinical situation.

What to do:


When you encounter an unexpected result, please repeat the test with a fresh qLabs® PT-INR test strip. If the same result is obtained, immediately consult with your local distributor.

NOTE:


- Please read the qLabs® PT-INR Test Strips Package Insert for more detail on the above content.

9.9 Finishing the Test

Remove and discard the used test strip into a sharps or other puncture proof container. It is recommended to power off the qLabs® ElectroMeter when it is not in use.

Before turning off the meter, press the “↩” button to return to the test mode.

Remove the test strip.


NOTE:

- If performing another test with a fresh test strip, please use a different finger for the new test.


10. Reviewing Previous Results

10.1 Reviewing the Previous Results


Follow these steps to review previous results that are stored in the meter memory.

In the test mode where the display prompts the user to insert a test strip.

Press the “MEM” button. The qLabs® ElectroMeter will enter the “Review/Delete” interface.

Select “Review” to enter Memory mode for reviewing previous results.

If no previous results exist, the meter will prompt “NO RECORD!”.

Whenever a test is successfully performed, the result is automatically stored in the meter memory along with the date of the test. The meter memory will store up to 200 results.

Entering Memory mode allows the user to:

01/01	PID0001
01/02	PID0002
02/01	PID0013
02/02	PID0015
	OK

Press the “↑” button to select previous results. Each page lists four results; use the functional the “↓” button to scroll down to review the next four results.

Press the “↶” button to return to the Test mode.

11.6 PT 0.8 INR	
8.2 QC	

Press the “OK” button to view the detailed information of the selected test including PT, INR results, and date. Use the “+” and “-” buttons to scroll through the detailed results. Press the “↶” button to exit.

CODE: 1505320	
LOT: 3052E0503	
01/05/2014 14:30	

Press the “↶” button to exit. Press the “→” button to view more detailed information and operation options.

10.1.1 Printing Patient Test Results

Press the “ ” button to enter “print setting” interface.


In print setting interface, press the “◀” and “▶” to select the number of test results to be printed. After selecting, press “OK” to confirm printing.

NOTE:


- A qLabs® eStation need to be connected to the qLabs® electroMeter to print selected results.

10.1.2 Deleting Patient Test Results

Press the “” button to enter “delete setting” interface.

Press the “OK” button to confirm or press “” button to exit.


If the option to delete the current result is selected, the “DELETED!” symbol will flash on the display for a short time, the subsequent results will be displayed.

If the option to delete all results is selected, the display will show “NO RECORD!”.


10.1.3 Uploading Patient Test Results


Press the “” button to upload the selected results to the qLabs® Data Manager.

10.2 Fast Deleting the Previous Results


Follow these steps to fast delete the previous results that are stored in the meter memory.


In the test mode where the display prompts the user to insert a test strip.

Press the “MEM” button. The qLabs® ElectroMeter will enter the “Review/Delete” interface.


Select “DELETE”. The qLabs® ElectroMeter will enter Fast Deleting Mode.

Press the “↓” button to scroll through the previous results, press “√” button to select the results which need to be deleted. Press “Delete” to delete the selected results. Press “Return” button to return to the previous interface.

When you pressed “DELETE” button, the meter will prompt you to confirm. Press “OK” to continue the deleting, press “Return” button to re-select results.

Once the deleting is finished, the “DELETED” symbol will flash on the display for a short time, and then return to the Fast Deleting Mode.

11. Quality Control

11.1 Onboard QC Test

The qLabs® ElectroMeter also measures the clotting time in the Quality Control (QC) zone. If the QC result fails to fall within a predetermined range, the qLabs® ElectroMeter will display an error code rather than give a possible erroneous PT-INR result.

This safety measure guards the user against situations in which the qLabs® PT-INR Test strip may have been subjected to very high temperatures or humidity (which could happen if the foil pouch is torn or punctured).

11.2 Liquid QC Test

Micropoint Biotechnologies, Inc. has optional liquid controls for qLabs® PT-INR system. The qLabs® PT-INR liquid controls are used with qLabs® PT-INR test strips to verify the qLabs® PT-INR system performance. It is recommended that two levels of liquid controls are used to perform the test.

The qLabs® PT-INR liquid controls can be purchased separately from Micropoint Biotechnologies, Inc. Please refer to the qLabs® PT-INR liquid controls package insert for detailed operation instructions.

Only use qualified qLabs® PT-INR Controls for QC test.

12. Maintenance

12.1 Care and Cleaning of Your Screen

Follow these steps for cleaning the qLabs® ElectroMeter:

NOTES:

No maintenance is required other than routine cleaning:

- Clean display screen with a clean damp cloth. If necessary, a mild detergent or disinfectant (such as 5% bleach solution or 70% isopropyl alcohol) may be used.
- Clean the area around the test strip guide with a swab or a pad that has been dampened with alcohol or 5% bleach solution.

12.2 Installing/Replacing Batteries

Follow these steps to install/replace the batteries:

1. Turn the qLabs® ElectroMeter upside down.
2. Remove the battery cover screw.
3. Remove the battery cover by pressing down on the battery cover release, remove the old batteries.
4. Replace with 4 standard AA batteries. Follow the direction for battery placement (+ and – ends) inside the battery compartment.
5. Fasten the battery cover. Tighten the battery cover screw.

NOTES:

- Turning the qLabs® ElectroMeter off manually when it is not being used will help conserve the battery power.
- Settings for date, time, past test results and reference range will be saved if batteries are replaced.
- If storing the ElectroMeter for an extended period, completely remove the batteries before the next use.

Dispose of the batteries in compliance with local regulations.

12.3 Servicing

All servicing and adjustment must be performed by Micropoint Biotechnologies, Inc.'s service engineers only.

13. Troubleshooting

ERROR CODE	DISCRIPTION	CORRECTIVE ACTIONS
	qLabs® ElectroMeter is not connected to qLabs® eStation.	Connect the qLabs® ElectroMeter to the qLabs® eStation properly.
	The qLabs® eStation is Out of thermal paper.	Replace the thermal paper.
NO COAG	The meter did not detect coagulation.	Retest with a new strip. If this error message appears again, the result must be checked using another method.
E001	Low battery power.	Replace batteries or use power adapter as a power source.
E002	Heating process timeout.	Turn off the meter and let it equilibrate at ambient temperature (10-35°C) for at least 5 minutes before performing additional tests.
E003	The ambient temperature is outside of the operational temperature range.	Move the meter and strip to a location that is within the operational range of 10C to 35C and repeat the test.
E004	Insufficient blood was added to the sample well.	Retest with a new strip and make sure a minimum of 10 micro liters of blood is correctly added to the sample well.
E005	Internal QC check out of specifications.	1. Make sure the strip is within the expiry date. 2. Make sure the strip code is entered correctly. 3. Retest with a new strip.
E006	External liquid QC check out of specifications.	1. Make sure the strip and liquid controls are within the expiry date. 2. Make sure the strip code is consistent with the code printed on the liquid controls vial label. 3. Retest with a new strip.

ERROR CODE	DISCRIPTION	CORRECTIVE ACTIONS
E007.x	PT calculation error.	Retest with a new strip.
E009	The INR is too low.	1. Make sure the strip code is entered correctly. 2. Retest with a new strip.
E010	The operation has timed out.	Please ensure the sample is added within 2 minute after the 'add sample' indicator is flashing
E011	You inserted a used test strip, or you applied a blood sample to the test strip before the meter warm-up was complete.	Retest with a new strip. Make sure to add sample within 2 minutes after the 'add blood' indicator is flashing.
E013	The strip is defective or im-properly operated when adding sample.	Retest using a new strip. Make sure the operation of adding sample is exactly correct according to the user manual guidance.
E014	The value of HCT may be out of range 30%-55%.	Patients with abnormal HCT value may be recommended to test the INR using another method.
E015	Bluetooth initialization error.	Turn off the meter and restart the meter. If the problem still exists, please contact the technician for further assistant.

14. Symbols

SYMBOLS	EXPLANATION	SYMBOLS	EXPLANATION
	In Vitro Diagnostics		Expiry Date
	Caution. Read Carefully		Fragile
	Avoid Liquid Spills		Biological Risk
	Separate Collection		DO Not Reuse
	Temperature Limitation		Manufacturer
	Authorized Representative in the European Community		Consult Instructions For Use
REF	Catalog Number		Serial Number
	CE Mark		Lot Number

15. Performance Characteristics & Product Specifications

15.1 Operating Condition

Temperature	between 10°C and 35°C (50°F and 95°F)
Humidity	10% to 90%
Atmospheric Pressure	700 hPa and 1060 hPa
Measurable Range	INR between 0.5 to 7.5

15.2 Product Specifications

CPU	32-Bit V2 ARM Cortex-M3
Display screen	128×96 LCD
Test memory	200 test results
Power	Battery: 4 X AA batteries Supply: Input: 100-240 VAC, 50-60 Hz Output: 7 VDC Power: 6 W
Size	135mm x 65mm x 34.5mm (L× W × H)
Weight	136 g (Without Batteries)
Warranty period	2 years

16. Warranty

Use of the qLabs® ElectroMeter

The qLabs® ElectroMeter (the “Meter”) is designed for use in monitoring patients on oral anticoagulant therapy. Proper adherence to the instructions in the User’s Manual and package insert is critical to proper operation.

WARNING: Failure to comply with the user manual could lead to inaccurate results and incorrect medication dosing which could result in injury or death.

Limited Warranty

Micropoint Biotechnologies, Inc. guarantees to the original purchaser of the Meter that the Meter is free from material defects in material and workmanship for two years from the date of purchase. This warranty does not guarantee the uninterrupted operation of the meter.

Micropoint Biotechnologies, Inc.’s only liability and purchaser’s only remedy under this warranty is that during the warranty period Micropoint Biotechnologies, Inc. shall replace or repair, at no charge, any Meter component with defects in material or workmanship. MICROPOINT BIOTECHNOLOGIES, INC. MAKES NO OTHER WARRANTIES AND EXPRESSLY EXCLUDES ANY IMPLIED WARRANTY OF MERCHANTABILITY, NONINFRINGEMENT OR FITNESS FOR A PARTICULAR USE.

The only user serviceable parts of the Meter are the battery, the cleaning door and the battery door. Tampering with any other portion of the Meter, abusing the Meter or using the Meter in a manner inconsistent with its user manual will void this warranty. This warranty does not apply to any component that is damaged by improper storage or accident or that is subject to alteration, misuse, tampering or abuse. Before returning any defective components, you must obtain a “Return Material Authorization” number and return instructions from Micropoint Biotechnologies, Inc. Technical Support by calling +86 755 86296766.

MICROPOINT BIOTECHNOLOGIES, INC. ENTIRE LIABILITY IN CONNECTION WITH THE METER, REGARDLESS OF THE LEGAL OR EQUITABLE BASIS OF ANY CLAIM, IS LIMITED TO THE PURCHASE PRICE OF THE METER. IN NO EVENT WILL MICROPOINT BIOTECHNOLOGIES, INC. BE LIABLE FOR ANY INCIDENTAL, INDIRECT, SPECIAL, CONSEQUENTIAL, OR PUNITIVE LOSS OR DAMAGES IN ANY WAY RELATED TO THE METER, OR ANY THIRD PARTY. CLAIMS, EVEN IF MICROPOINT BIOTECHNOLOGIES, INC. OR ITS DISTRIBUTORS ARE AWARE OF THE POSSIBILITY OF SUCH CLAIMS OR DAMAGES AND NOTWITHSTANDING THE ESSENTIAL PURPOSE OF ANY REMEDY.

17. Contact Information

Customer Service and Technical Support

Available 8:30 AM to 5:30 PM (Beijing Standard Time), Monday-Friday.

Tel: +86 755 86296766

Micropoint Biotechnologies, Inc.

6F, No.3, Industry 5 Road, Shekou,
Shenzhen, Guangdong 518067, China

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©2015 Micropoint Biotechnologies, Inc. All rights reserved. Printed in China.
qLabs® and Micropoint® are registered trademarks of Micropoint Bioscience, Inc.